



**TZ-28**  
 Titanium Rotor  
 1,350 ml  
 Zonal  
 Maximum Speed: 20,000 rpm  
 K Factor: 612  
 Rotor Mass: 7.80 kg

	RCF (g)	R (cm)
Max.	42,580	9.53
Min.	16,174	3.62

		RC-28S	RC-26 Plus	RC-24 <sup>1</sup>	RC-5C Plus RC-5C <sup>1</sup>	RC-5B Plus RC-5B <sup>1</sup> RC-5 <sup>1</sup>	RC-2B <sup>1</sup> RC-2 <sup>1</sup>	SS-3 <sup>1</sup> SS-4 <sup>1</sup>	RC-3B <sup>1</sup> RC-3B Plus RC-3C <sup>1</sup> RC-3C Plus
<b>TZ-28</b>	<b>Max Speed (rpm)</b>	20,000	20,000	20,000	20,000	20,000	—	—	—
	<b>Max RCF (g)</b>	42,580	42,580	42,580	42,580	42,580	—	—	—

<b>Catalog No. 52358</b>	<b>Required for Zonal Operation</b>	<b>Required for Sealed Operation</b>
Includes: TZ-28 Rotor Core Septa <b>49609</b> O-ring for Core Septa 49609 (12 supplied) <b>66193</b> O-ring for Rotor Cover (2 supplied) <b>65755</b> Rotor Cover Wrench <b>49724</b> T-wrench <b>49026</b> Rotor Holding Fixture <b>52540</b> Sample Inject Needle <b>66026</b> Lubricant <b>61556</b> Feed Adapter <b>49007</b> O-ring for Feed Adapter <b>63037</b>	Drive Adapter Conversion Kit <b>12344</b> Includes Manual <b>49618</b> Zonal Operation Kit <b>12343</b> RC-28S Tubing Support Bracket (Only necessary for update of pre-1989 zonal kit) <b>59637</b>	Drive Adapter Conversion Kit <b>12344</b> Zonal Operation Kit <b>12343</b> Sealing Assembly <b>12273</b> RC-28S Tubing Support Bracket (Only necessary for update of pre-1989 zonal kit) <b>59637</b>

<sup>(1)</sup>Instrument no longer available.