



TZ-28/GK System

Titanium
 1,350 ml
 Continuous Flow
 Max. Collectible Precipitate: 800 ml
 Max. Flow Rate: 600 ml/min.
 Maximum Speed: 19,000 rpm
 K Factor: 678
 Rotor Mass: 7.80 kg

	RCF (g)	R (cm)
Max.	38,428	9.53
Min.	14,597	3.62

		RC-28S	RC-26 Plus	RC-24 ¹	RC-5C Plus RC-5C ¹	RC-5B Plus RC-5B ¹ RC-5 ¹	RC-2B ¹ RC-2 ¹	SS-3 ¹ SS-4 ¹	RC-3B ¹ RC-3B Plus RC-3C ¹ RC-3C Plus
TZ-28/GK	Max Speed (rpm)	19,000	19,000	19,000	19,000	19,000	n/a	n/a	n/a
	Max RCF (g)	38,428	38,428	38,428	38,428	38,428	n/a	n/a	n/a

B-2

Catalog No. 52358	Required Accessories	Optional Accessories
Includes: TZ-28 Rotor Core Septa 49609 O-ring for Core Septa 49609 (12 supplied) 66193 O-ring for Rotor Cover (2 supplied) 65755 Rotor Cover Wrench 49724 T-wrench 49026 Rotor Holding Fixture 52540 Sample Inject Needle 66026 Lubricant 61556 Feed Adapter 49007 O-ring for Feed Adapter 63037	Drive Adapter Conversion Kit (Not Necessary if TZ-28 Zonal SUPERSPEED Already Purchased) 12344 Includes Manual 49618 Continuous Flow Distributor Adapter for Inlet-Outlet Assembly 49633 Inlet-Outlet Assembly for Distributor Adapter 49089	Bearings Kit, Replacement 12779 Syringe Assembly for Bearing Lubrication 13127 High Speed Lift Assembly 13086 RC-28S Tubing Support Bracket (Updates Kits for use in RC-28S; included with Inlet-Outlet Assembly 49089) 59633

⁽¹⁾Instrument no longer available.