



SS-34/KSB System

Aluminum Rotor
 8-Place, 50 ml
 Fixed-Angle: 34°, Continuous Flow
 Max. Collectible Precipitate: 280 ml
 Max. Flow Rate: 600 ml/min.
 Maximum Speed: 20,000 rpm
 K Factor: 750
 Rotor Mass: 6.8 kg
 Autoclavable: 121°C

	RCF (g)	R (cm)
Max.	47,808	10.70
Min.	14,610	3.27

		RC-28S	RC-26 Plus	RC-24 ¹	RC-5C Plus RC-5C ¹	RC-5B Plus RC-5B ¹ RC-5 ¹	RC-2B ¹ RC-2 ¹	SS-3 ¹ SS-4 ¹	RC-3B ¹ RC-3B Plus RC-3C ¹ RC-3C Plus
SS-34/KSB	Max Speed (rpm)	20,000	20,000	20,000	20,000	20,000	20,000	n/a	n/a
	Max RCF (g)	47,808	47,808	47,808	47,808	47,808	47,808	n/a	n/a

Catalog No. 14005	Optional Accessories
Includes: Complete SS-34/KSB System Manual 14008 <u>Note: SS-34 Rotor 28020 not included.</u>	High Speed Lift Assembly: Designed to reduce or eliminate the excess foaming sometimes encountered when operating at high speeds, particularly with certain nutrient broths or ammonium sulfate solutions. 13086 Four-Channel Distributor for KSB System. Allows collection in four tubes. 13109 Bearings Kit for KSB Continuous Flow System, Replacement. 12779 Tubing Support Bracket for RC-2, RC-2B, RC-5, RC-5B, RC-5C; Converts KSB-3 to KSB-R 14000 Rotor-Locking Stud 13017 Tubing Support Bracket for RC-28S 59633 (Updates previously purchased Kits for use in RC-28S; now included with KSB Kit 14005 .)

⁽¹⁾Instrument no longer available.

B-2