



SLC-250T
SUPER-LITE® Composite Rotor
 6-Place, 250 ml
 Tube Dimensions: 61 × 153 mm
 Fixed-Angle: 23°
 Maximum Speed: 11,500 rpm
 K Factor: 2,484
 Rotor Mass: 7.02 kg
 Autoclavable: 121°C

	RCF (g)	R (cm)
Max.	20,238	13.70
Min.	5,525	3.74

Used in SUPER T 21 centrifuge.

Catalog No. 76902		Accessories
Includes: SLC-250T Rotor Vacuum Grease 65937 Rotor Cover Assembly 76904	Rotor Cover O-Ring (6 supplied) 92524 Rotor Stand 76911 Instruction Manual 76905	Rotor Extractor 76818 ¹ Retaining Ring Pliers 65614 ¹ Software Update Kit 76026

TUBES AND BOTTLES FOR SLC-250T

CAPACITY (ml)		DESCRIPTION	MAX SPEED (rpm)	DIMENSIONS (mm)	TUBE CAT. NO. (Qty/Pkg)	SEALING ASSEMBLY*		ADAPTER REQ'D		NOTES
NOMINAL	FILL					CAT. NO. (Qty/Pkg)	DESCRIPTION	CAT. NO. (Qty/Pkg)	NO. OF PLACES	
250	270	Polycarbonate DRY-SPIN Bottle	11,500	61 × 146	03098 (12)	Included with Bottle	Polypropylene, Sealing	—	—	Replacement O-Ring sealing assembly 03545 (24/pkg). Replacement Caps 06841 (12/pkg).
250	280	Polypropylene DRY-SPIN Bottle	11,500	61 × 146	03069 (12)	Included with Bottle	Polypropylene, Sealing	—	—	Replacement O-Ring sealing assembly 03545 (24/pkg). Replacement Caps 06841 (12/pkg).
250	200	Polyallomer Oak Ridge Bottle	11,500	60 × 120	54477 (6)	54421 (1)	Delrin® Screw Cap	12002 (1)	1	
250	225	Polypropylene Oak Ridge Bottle	11,500	61 × 124	03937 (12)	03278 (12)	Polypropylene, Sealing	—	—	Nonsealing Cap 03927 (12/pkg) available.
250	240	Polycarbonate Oak Ridge Bottle	11,500	61 × 124	03939 (12)	03278 (12)	Polypropylene, Sealing	—	—	Nonsealing Cap 03927 (12/pkg) available.
150	144	Polycarbonate Thick-Walled Tube	11,500	45 × 132	03143 (25)	—	—	00458 (1)	1	
150	120	Corex®	5,000	53 × 132	03229 (1)	Included with Bottle	Polypropylene, Nonsealing	00372 (1)	1	Bottle no longer available
125	100	Corex®	5,000	53 × 118	03231 (1)	Included with Bottle	Polypropylene, Nonsealing	00372 (1)	1	Bottle no longer available
100	100	Pyrex® Bottle	3,000	44 × 137	03230 (1)	Included with Bottle	Polypropylene, Nonsealing	00371 (1)	1	
50	—	Tissue Culture Tube, Conical	—	30 × 115	—	—	—	03072 (1)	1	
35	—	Polyallomer Nalgene® Conical Oak Ridge Tube	—	29 × 105	—	—	—	03072 (1)	1	
35	—	Polycarbonate Nalgene® Conical Oak Ridge Tube	—	29 × 105	—	—	—	03072 (1)	1	

B-3

⁽¹⁾These tools are included with the SUPER T 21 instrument.

(continued on next page)

*Sealing Assembly replacement parts may be available. See section E for availability.

All trademarks are the property of their respective owners.

