

SORVALL® LOWSPEED AND TABLETOP ROTOR TO CENTRIFUGE CROSS-REFERENCE

CENTRIFUGES										
MAXIMUM SPEED (rpm)										
LOWSPEED AND TABLETOP ROTORS	FLOOR MODEL		TABLETOP							
	RC-3BP RC-3B PLUS RC-3B ¹ RC-3C PLUS RC-3C ¹	RC-3 ¹	RT-6000D ¹ RT-6000B ¹ T-6000B ¹ T-6000D	RT-6000 ¹ T-6000 ¹ GLC-4 ¹	GLC-2B ¹	GLC-3 ¹	TC-6	RMC-14	MC-12 MC-12V MC-12C ¹	RT-7
F-12/M.18									12,000	
FA-MICRO								14,000	12,000	
M		5,250			6,000	3,400				
A-384			6,000	6,000						
SP/X		5,250			6,000	3,400				
A-500			6,000	6,000						
A/S-400			4,380							
LA/S-400	6,000									
SL-50RT										7,000
HL-4					3,200	3,400				
H-400							3,500			
H-1000B			3,200	3,200						
RTH-250										3,200
RTH-750										4,000
H-2000B	2,800									
HBB-12	2,900									
HG-4L ³	5,250 ²	5,250								
H-4000 ³	5,250									
H-6000A	5,000									
HL-2B	3,300									
HBB-6	5,000									

⁽¹⁾Centrifuge no longer available.

⁽²⁾Requires Adapter Conversion Kit Cat. No. 12300.

⁽³⁾Rotor no longer available.