



6 x 94 ml Rotor
 Aluminum
 6- Place, 94 ml
 Tube Dimensions: 38.1 x 102 mm
 Fixed Angle: 26°
 Maximum Speed: 12,000 rpm
 K-Factor: 829
 Rotor Weight (Kg/lb): 5.5 / 12.1
 Autoclaveable: 121°C

	R(cm)	RCF (g)
Max	10.1	16,260
Min	6.3	10,142

For use in the Legend T/RT and the Multifuge 3/3R
 (Rotor can also be used in the Biofuge / Contifuge Stratos
 @ 15,000 rpm / 25,410 x g)

6 x 94 ml Rotor: 75003334

INCLUDES:

Rotor Cover O-Ring: 20290628
 Trunnion Grease: 70006692 (Supplied with centrifuge)
 O-ring Grease: 76003500 (Supplied with centrifuge)

CAPACITY (ml)		Description	Max Speed (rpm)	Dimensions (mm)	Tube Cat. No. (Qty/Pkg)	SEALING ASSEMBLY		ADAPTER REQUIRED		NOTES
Nominal	Fill					Cat. No. (Qty/Pkg)	Description	Cat. No. (Qty/Pkg)	No. of places	
94	70	Polycarbonate tube	12,000	38.1 x 101.6	76002810 (1)	75001568 (1)	Aluminum with O-ring	--	--	Sealing caps included with tubes
94	70	Polyallomer tube	12,000	38.1 x 101.6	76002872 (1)	75001568 (1)	Aluminum with O-ring	--	--	Sealing caps included with tubes
94	94	Polycarbonate Nalgene® UltraBottles	12,000	38 x 102	3432-0070*	3430 - 3870*	Polyphenylene Oxide	--	--	Autoclaveable
						3431 - 3870*	Aluminum			
85	85	Polycarbonate Oakridge Bottle	12,000	38 x 105	75002976 (10)	included		--	--	
50	46	Polypropylene Flanged Tube	12,000	29 x 102	03147 (25)	03268 (25)	Polypropylene Snap-On	75003014	1	
50	48	Polycarbonate Flange Tube	12,000	29 x 102	03146 (25)	03268 (25)	Polypropylene Snap-On	75003014	1	
50	56	Polyallomer Thin-Walled Tube	12,000	28 x 104	03139 (25)	--	--	75003014	1	
50	41	Polypropylene Oakridge Bottle	12,000	28 x 103	03530 (25)	03535 (25)	Polypropylene Sealing	75003014	1	
50	38	Polycarbonate Oakridge Bottle	12,000	28 x 103	03480 (25)	03535 (25)	Polypropylene Sealing	75003014	1	
50	40	Polycarbonate Screw Cap	12,000	28.5 x 105	75002963 (10)	75002967 (2) incl.	Polypropylene Sealing	75003014	1	
50	40	Polypropylene Copolymer Screw Cap	12,000	28.5 x 105	75002964 (10)	75002967 (2) incl.	Polypropylene Sealing	75003014	1	
50	45	Teflon Screw Cap	12,000	28.5 x 105	75002966 (2)	75002967 (2) incl.	Polypropylene Sealing	75003014	1	
50	50	Stainless Steel Tube	12,000	28 x 101	00579 (1)	00518 (1)	Stainless Steel Sealing	75003014	1	
50	50	Stainless Steel Flange Tube	12,000	28 x 101	00517 (1)	00518 (1)	Stainless Steel Sealing	75003014	1	
38		Polypropylene Copolymer oakridge bottle	12,000	25 x 104	75003086			75003094	1	
35	--	Polyallomer Nalgene® Conical Oak Ridge Tube	--	29 x 105	--	--	--	75003096	1	

* Tube must be ordered from Nalgene® (1-800-625-4327)

CAPACITY (ml)		Description	Max Speed (rpm)	Dimensions (mm)	Tube Cat. No. (Qty/Pkg)	SEALING ASSEMBLY		ADAPTER REQUIRED		NOTES
Nominal	Fill					Cat. No. (Qty/Pkg)	Description	Cat. No. (Qty/Pkg)	No. of places	
35	--	Polycarbonate Nalgene® Conical Oak Ridge Tube	--	29 x 105	--	--	--	75003096	1	
30	--	Polycarbonate Oak Ridge Bottle	12,000	25 x 94	03719 (25)	03276 (25)	Polypropylene Sealing	75003094	1	
30	--	Polypropylene Oak Ridge Bottle	12,000	25 x 94	03718 (25)	03276 (25)	Polypropylene Sealing	75003094	1	
30	31	Corex® Bottle	12,000	24 x 105	03163 (10)	--	--	75003014 + 00368	1	
16	16	Polypropylene Flanged Tube	12,000	18 x 100	03244 (50)	03299 (50)	Polypropylene Non-sealing	76002906	1	
16	17	Polycarbonate Flange Tube	12,000	18 x 100	03243 (50)	03269 (50)	Polypropylene Non-sealing	76002906	1	
16		Polycarbonate thick walled tube	12,000	18 x 100	75002950 (10)	75002957 (10)	75002957 Polypropylene stoppers	76002906	1	
16		PP Copolymer thick walled tube	12,000	18 x 100	75002951	75002957 (10)	75002957 Polypropylene stoppers	76002906	1	
15	18	Corex® Bottle	12,000	18 x 102	03286 (25)	--	--	76002906	1	
12	--	Falcon snap cap tube (# 2059)	--	17 x 100	Please contact you local Kendro Laboratory Products representative for more information.					
12	10	Polypropylene Flanged Tube	12,000	16 x 100	03116 (50)	03266 (50)	Polypropylene Non-sealing	75003014 + 00402	1	Both adapters required
12	11	Polycarbonate Flange Tube	12,000	16 x 100	03115 (50)	3267 (50)	Polypropylene Non-sealing	75003014 + 00402	1	
10	11	Polyallomer Thin-Walled Flange Tube	12,000	13 x 100	03124 (50)	--	--	75003092	2	
10	--	Glass Tube	--	13 x 100	--	--	--	75003092	2	
10	9.1	Polypropylene Oakridge Bottle	12,000	16 x 83	03929 (50)	03279 (25)	Polypropylene Sealing	75003093	2	
10	9.7	Polycarbonate Oakridge Bottle	12,000	16 x 83	03020 (25)	03279 (25)	Polypropylene Sealing	75003093	2	
7	7	Polypropylene Flanged Tube	12,000	13 x 100	03121 (50)	03265 (50)	Polypropylene Non-sealing	75003092	2	
7	7	Polycarbonate Flange Tube	12,000	13 x 100	03120 (50)	03265 (50)	Polypropylene Non-sealing	75003092	2	
5	--	Falcon snap cap tube (# 2063)	--	12 x 75	Please contact you local Kendro Laboratory Products representative for more information.					
4	4	Polypropylene Flanged Tube	12,000	11 x 75	03105 (50)	03264 (50)	Polypropylene Non-sealing	75003091	4	
4	4	Polycarbonate Flange Tube	12,000	11 x 75	03104 (50)	03264 (50)	Polypropylene Non-sealing	75003091	4	
1.5	--	Microtubes	12,000	--	--	--	--	76002905	4	