



HIGHconic® Rotor

Aluminum (?)
 6- Place, 50 ml
 Tube Dimensions: 30 x 117 mm
 Fixed Angle: 45°
 Maximum Speed: 8,500 rpm
 K-Factor: 2,542
 Rotor Weight (Kg/lb): 2.8 / 6.2
 Autoclaveable: 121°C

	R(cm)	RCF (g)
Max	12.4	10,016
Min	6.0	4,847

**For use in the Legend T/RT / Multifuge 3/3R /
 Biofuge Stratos / Contifuge Stratos**

HIGHconic Rotor: 75003046

INCLUDES:

- 50-ml conical adapter: 75003096
- Rotor Cover O-Ring: 20290628
- Trunnion Grease: 70006692 (Supplied with centrifuge)
- O-ring Grease: 76003500 (Supplied with centrifuge)

CAPACITY (ml)		Description	Max Speed (rpm)	Dimensions (mm)	Tube Cat. No. (Qty/Pkg)	SEALING ASSEMBLY		ADAPTER REQUIRED		NOTES
Nominal	Fill					Cat. No. (Qty/Pkg)	Description	Cat. No. (Qty/Pkg)	No. of places	
50	46	Polypropylene Flanged Tube	8,500	29 x 102	03147 (25)	03268 (25)	Polypropylene Snap-On	75003014	1	
50	48	Polycarbonate Flange Tube	8,500	29 x 102	03146 (25)	03268 (25)	Polypropylene Snap-On	75003014	1	
50	56	Polyallomer Thin-Walled Tube	8,500	28 x 104	03139 (25)	--	--	75003014	1	
50	41	Polypropylene Oakridge Bottle	8,500	28 x 103	03530 (25)	03535 (25)	Polypropylene Sealing	75003014	1	
50	38	Polycarbonate Oakridge Bottle	8,500	28 x 103	03480 (25)	03535 (25)	Polypropylene Sealing	75003014	1	
50	40	Polycarbonate Screw Cap	8,500	28.5 x 105	75002963 (10)	75002967 (2) incl.	Polypropylene Sealing	75003014	1	
50	40	Polypropylene Copolymer Screw Cap	8,500	28.5 x 105	75002964 (10)	75002967 (2) incl.	Polypropylene Sealing	75003014	1	
50	45	Teflon Screw Cap	8,500	28.5 x 105	75002966 (2)	75002967 (2) incl.	Polypropylene Sealing	75003014	1	
50	50	Stainless Steel Tube	8,500	28 x 101	00579 (1)	00518 (1)	Stainless Steel Sealing	75003014	1	
50	50	Stainless Steel Flange Tube	8,500	28 x 101	00517 (1)	00518 (1)	Stainless Steel Sealing	75003014	1	
50	--	Conical Tissue Culture Tube	--	30 x 115	--	--	--	75003096	1	Adapter supplied with rotor
38	--	Polypropylene Copolymer oakridge bottle	8,500	25 x 104	75003086			75003094	1	
35	--	Polyallomer Nalgene® Conical Oak Ridge Tube	--	29 x 105	--	--	--	75003096	1	Adapter supplied with rotor
35	--	Polycarbonate Nalgene® Conical Oak Ridge Tube	--	29 x 105	--	--	--	75003096	1	
30	--	Polycarbonate Oak Ridge Bottle	8,500	25 x 94	03719 (25)	03276 (25)	Polypropylene Sealing	75003094	1	

CAPACITY (ml)		Description	Max Speed (rpm)	Dimensions (mm)	Tube Cat. No. (Qty/Pkg)	SEALING ASSEMBLY		ADAPTER REQUIRED		NOTES
Nominal	Fill					Cat. No. (Qty/Pkg)	Description	Cat. No. (Qty/Pkg)	No. of places	
30	--	Polypropylene Oak Ridge Bottle	8,500	25 x 94	03718 (25)	03276 (25)	Polypropylene Sealing	75003094	1	
30	31	Corex® Bottle	8,500	24 x 105	03163 (10)	--	--	75003014 + 00368	1	
16	16	Polypropylene Flanged Tube	8,500	18 x 100	03244 (50)	03299 (50)	Polypropylene Non-sealing	76002906	1	
16	17	Polycarbonate Flange Tube	8,500	18 x 100	03243 (50)	03269 (50)	Polypropylene Non-sealing	76002906	1	
16		Polycarbonate thick walled tube	8,500	18 x 100	75002950 (10)	75002957 (10)	75002957 Polypropylene stoppers	76002906	1	
16		PP Copolymer thick walled tube	8,500	18 x 100	75002951	75002957 (10)	75002957 Polypropylene stoppers	76002906	1	
15	--	Conical Tissue Culture Tube	--	17 x 120	--	--	--	75003095	1	
15	18	Corex® Bottle	8,500	18 x 102	03286 (25)	--	--	76002906	1	
12	--	Falcon snap cap tube (# 2059)	--	17 x 100	Will require a custom adapter. Please contact you local Kendro Laboratory Products representative for more information.					
12	10	Polypropylene Flanged Tube	8,500	16 x 100	03116 (50)	03266 (50)	Polypropylene Non-sealing	75003014 + 00402	1	Both adapters required
12	11	Polycarbonate Flange Tube	8,500	16 x 100	03115 (50)	3267 (50)	Polypropylene Non-sealing	75003014 + 00402	1	
10	11	Polyallomer Thin-Walled Flange Tube	8,500	13 x 100	03124 (50)	--	--	75003092	2	
10	--	Glass Tube	--	13 x 100	--	--	--	75003092	2	
10	9.1	Polypropylene Oakridge Bottle	8,500	16 x 83	03929 (50)	03279 (25)	Polypropylene Sealing	75003093	2	
10	9.7	Polycarbonate Oakridge Bottle	8,500	16 x 83	03020 (25)	03279 (25)	Polypropylene Sealing	75003093	2	
7	7	Polypropylene Flanged Tube	8,500	13 x 100	03121 (50)	03265 (50)	Polypropylene Non-sealing	75003092	2	
7	7	Polycarbonate Flange Tube	8,500	13 x 100	03120 (50)	03265 (50)	Polypropylene Non-sealing	75003092	2	
5	--	Falcon snap cap tube (# 2063)	--	12 x 75	Will require a custom adapter. Please contact you local Kendro Laboratory Products representative for more information.					
4	4	Polypropylene Flanged Tube	8,500	11 x 75	03105 (50)	03264 (50)	Polypropylene Non-sealing	75003091	4	
4	4	Polycarbonate Flange Tube	8,500	11 x 75	03104 (50)	03264 (50)	Polypropylene Non-sealing	75003091	4	
1.5	--	Microtubes	8,500	--	--	--	--	76002905	4	