



F-12/M.18
 Aluminum Rotor
 18-Place, 1.5 ml
 Tube Dimensions: 11 × 38 mm
 Fixed-Angle: 45°
 Maximum Speed: 12,000 rpm
 K Factor: 732
 Rotor Mass: 0.81 kg

	RCF (g)	R (cm)
Max.	12,128	7.54
Min.	7,994	4.97

For use in MC-12V, MC-12C¹, and MC-12 centrifuges

Catalog No. 73002	Accessories
Includes: F-12/M.18 Rotor Rotor Cover 73164 Rotor Tie Down Stud 73166 Washer for Rotor Stud 90812 Rotor Instruction Manual 73006 Rotor Wrench 73031	Rotor Extractor Tool 12935

TUBES

NOMINAL CAPACITY (ml)	DESCRIPTION	ADAPTER REQ'D.		MAX RADIUS	MAX g FORCE
		CAT. NO. (Qty/Pkg)	NO. OF PLACES		
1.5	Microtube	—	—	7.22	11,619
0.5	Microtube	03722 (18)	1	6.18	9,958
0.4	Microtube	03724 (18)	1	7.54	12,127
0.25	Microtube	03724 (18)	1	6.13	9,864

⁽¹⁾Instrument no longer available.