



HBB-6
 Stainless Steel Rotor Body
 Aluminum Swinging Buckets
 6-Place, 1400 ml
 Cavity Dimensions:
 90W × 133L × 140D mm
 Maximum Speed: 5,000 rpm
 Rotor Mass: 25.55 kg

	RCF (g)	R (cm)
Max.	7,129	25.53

For use in RC3BP®, RC-3C Plus, RC-3C,¹ RC-3B Plus and RC-3B¹ centrifuges

Catalog No. 11753	Accessories
Includes: HBB-6 Rotor Rotor Cover Assembly 11237 (6) Rotor Buckets, Anodized Aluminum, 11758 Rotor Locking Tool 11246 Rotor Cleaning Kit 12259 Balancing Discs, Rubber (1 pkg, 6/pkg) 00335 Instruction Manual 11243	Microplate Carriers (set of 2) 11267 Rotor Buckets 11754 (2/pkg) Bucket Liners, divided cavity (2/pkg) 11751 Bucket Liners, single cavity (2/pkg) 11756 Rotor Extractor 11445 Conversion Kit 11758 Adapter (1/pkg) 01094 Adapter (1/pkg) 01098
 00335	 11267
	 01098
	 01094
	 11751
	 11756

BLOOD BAGS

CAPACITY (ML)		DESCRIPTION	MAX SPEED (RPM)	DIMENSIONS (MM)	TUBE CAT. NO. (QTY/PKG)	SEALING ASSEMBLY		ADAPTER REQ'D		NOTES
NOMINAL	FILL					CAT. NO. (QTY/PKG)	DESCRIPTION	CAT. NO. (QTY/PKG)	NO. OF BAGS/PLACES	
—	—	Platelet Rich Plasma or Buffy Coat	—	—	—	—	—	11756 ² (2)	2-4	Depends on volume.
—	400-450	Blood Bag, Single and Double ³	—	—	—	—	—	11751 (2)	2	Divided liner.
—	400-450	Whole Blood Bag, Double or Triple	—	—	—	—	—	01094 (1)	1	
—	450-550	Whole Blood Bag, Quad or Quint	—	—	—	—	—	01098 (1)	1	
—	450-550	Whole Blood Bag, Quad ³	—	—	—	—	—	01094 (1)	1	

⁽¹⁾Instrument no longer available.
⁽²⁾Bucket liner is optional.
⁽³⁾Use when no Adsol™ is present in satellite bag.

All trademarks are the property of their respective owners.