





HBB-12
 Stainless Steel Rotor Body
 Aluminum Buckets
 6-Place
 2 Quad Blood Bags/Place or
 3 Satellite Bag Sets/Place
 Swinging-Bucket
 Maximum Speed: 2,900 rpm
 Rotor Mass: 26.7 kg

	RCF (g)	R (cm)
Max.	2,462	26.21

For use in **RC3BP®**, **RC-3C Plus**, **RC-3C**,¹ **RC-3B Plus** and **RC-3B**¹ centrifuges

Catalog No. 11276	Accessories
Includes: HBB-12 Rotor Rotor Buckets, Anodized Aluminum 11297 Rotor Cleaning Kit 12259 Balancing Discs, Rubber (1 pkg, 6/pkg) 00335 Instruction Manual 11283	Rotor Extractor 11445 Liner 11284 is a sleeve for the buckets, for ease of loading and unloading blood bags.
 00335	 11284

TUBES, BOTTLES

CAPACITY (ml)		DESCRIPTION	MAX SPEED (rpm)	DIMENSIONS (mm)	ADAPTER REQ'D	
NOMINAL	FILL				CAT. NO. (Qty/Pkg)	NO. OF BAGS/PLACES
—	450	Blood Bag Triple	—	—	11284 (1)	2

⁽¹⁾Instrument no longer available.