



HL-2B
 RIA Tube Rack Rotor
 Aluminum Rotor Body and Carriers
 2-Place, Multiple Rack
 Maximum Speed: 3,300 rpm
 Rotor Mass: 10.4 kg

	RCF (g)	R (cm)
Max.	2,024	16.64

For use in **RC3BP®**,
RC-3C Plus, **RC-3C**,¹
RC-3B Plus and **RC-3B**¹
 centrifuges

Catalog Numbers

- 11421:** HL-2B Tube Rack Rotor with Trays and Clips for Beckman System
- 11417:** HL-2B Tube Rack Rotor with Trays and Clips for LKB System
- 11410:** HL-2B Tube Rack Rotor with Trays and Clips for Micromedics System
- 11418:** HL-2B Tube Rack Rotor with Trays and Clips for Packard System
- 11411:**² HL-2B Tube Rack Rotor with Trays and Clips for Tracor/TM Analytic System
- 11445:** Rotor Extractor Kit, Optional

Replacement Parts

RACKS (NOT MANUFACTURED BY SORVALL®)	RACK CAPACITY (NO. OF TUBES)	NUMBER OF RACKS THE HL-2B ACCEPTS	TOTAL NUMBER OF TUBES CARRIED	CATALOG NUMBERS FOR HL-2B Rotor Parts ³	
				TRAY (2/pkg)	ADAPTER CLIPS (4/pkg)
Beckman	10	12	120	42578	42579
LKB	11	16	176	42557	42567
Micromedics	14 or 20	12	168 or 240	42550	42560
Packard	12	6	72	42558	42568
Tracor/TM Analytic	24	6	144	42551	42561
Tracor/TM Analytic	36	4	144	42551	42561

⁽¹⁾Instrument no longer available.

⁽²⁾Requires Support Block or rack: P/N 000003615. Purchase separately and submit to SORVALL for modification.

⁽³⁾Also for use in HL-2 rotor in RC-3 Centrifuge.