


ATC Probe 917004



- **ATC Temperature Sensor**
- **Durable stainless steel construction**
- **2.5 mm ATC plug**

917004

Lot Code:

Thermo
SCIENTIFIC

OPERATION INSTRUCTIONS

Thank you for selecting this ATC probe.
Temp: 0-100 °C. Dims (mm): 148 L x 3 dia.

DURING MEASUREMENT: To avoid contamination, rinse probe with deionized water before and between measurement. Stir sample to ensure homogeneous mixture.

WARRANTY: This probe is warranted for a period of 6 months from the date of purchase. If replacement is necessary within the 6 month period and the electrode issue is not the result of abuse or misuse, please contact your local distributor.

RETURN OF ITEM: Authorization must be obtained before returning items for any reason. Note: Prices and design are subject to change without notice.

Made in Singapore

Distributed by Thermo Fisher Scientific

Environmental Instruments **Water Analysis Instruments**

166 Cummings Center, Beverly, MA 01915 USA
Toll Free: 1-800-225-1480 • Dom. Fax: 1-978-232-6015
Tel: 1-978-232-6000 • Int'l Fax: 978-232-6031

www.thermo.com/water