



Partners in Protection

Trust Thermo Scientific laboratory equipment to protect you and your valuable cell cultures.



Thermo Scientific 1300 Series A2 Biological Safety Cabinets

Enhanced Safety

Dual DC motors continuously control and verify cabinet airflow, while the performance monitor clearly displays inflow/downflow values.

Exceptional Ergonomics

Features slanted sash for enhanced posture and reach, coated back walls to prevent glare and low noise level for better concentration.

Environmental Protection

Reduces carbon footprint up to 80% from competitive units due to the energy savings of dual DC motors. Visit www.thermo.com/1300 to read our new energy article and to see a virtual 3D demo of the 1300.

Thermo Scientific NAPCO® CO₂ Incubators

Secured

Choose with confidence: water jacket models for enhanced temperature stability, or direct heat units that offer a 140°C sterilization cycle to easily clean the chamber.

Protected

Patented HEPA airflow system creates a Class 100 environment after door closing, eliminating potential contamination risk.

Thriving

In-chamber sensors and gentle air circulation assure a stable environment for optimal cell growth. Models with O₂ control are available, along with IR sensors and humidity displays.




► **Purchase our 1300 Series A2 biological safety cabinet or our NAPCO CO₂ incubator and receive your choice of a FREE accessory.**

More details on the back.

Get an **extra bonus** when you protect with the best!

► Purchase a Thermo Scientific 1300 Series A2 biological safety cabinet or a Thermo Scientific NAPCO CO₂ incubator and **receive a FREE accessory!**



Buy our 1300 Series A2* Biological Safety Cabinet

Choose 1 of two **FREE accessories** —

Foot Rest > Value of \$210

Our 1300 Series A2 packages include an adjustable stand, arm rest and factory-installed UV light. Choose the foot rest to increase your comfort level while working in the cabinet.

Vacuum Valve > Value of \$103

Provides an access to the vacuum in your cabinet.

*Model Numbers 1385 or 1387



Buy our NAPCO CO₂ Incubator

Choose 1 of two **FREE accessories** —

2 HEPA Replacement Kits > Value of \$161

Each kit includes one HEPA filter and two smaller in-line filters.

Our exclusive HEPA filtration system provides Class 100 air quality in 5 minutes. Protects the cells from airborne contamination by filtering the chamber air every 60 seconds.

Gas Guard > Value of \$771

The gas guard is an external CO₂ tank switch to protect your CO₂ levels.



► Buy both products together and choose **an additional FREE accessory!**

ORDERING INFORMATION

Thermo Scientific 1300 Series A2 Biological Safety Cabinet

Cat No.	Description
13-261-304	1385 A2 4-foot Biosafety Cabinet Package with adjustable height stand, UV light, and arm rest
13-261-306	1387 A2 6-foot Biosafety Cabinet Package with adjustable height stand, UV light, and arm rest

Offers expire on 12/15/2008. Products must be purchased between 7/1/2008 -12/15/2008. Special offers must be redeemed by 12/15/2008. Valid in the U.S. only. Not valid with any other offers, contracts or agreements. Void where prohibited. All offers based on U.S. List prices.

Thermo Scientific NAPCO CO₂ Water Jacketed Incubator

Cat No.	Description
15-497-001	Incubator with TC sensor
15-497-002	Incubator with TC sensor, RH display
15-497-003	Incubator with IR sensor
15-497-004	Incubator with IR sensor, RH display
15-497-005	Incubator with TC sensor, O ₂ control
15-497-006	Incubator with TC sensor, RH display, O ₂ control
15-497-007	Incubator with IR sensor, O ₂ control
15-497-008	Incubator with IR sensor, RH display, O ₂ control

Thermo Scientific NAPCO CO₂ Direct Heat Incubator

Cat No.	Description
15-497-009	Incubator with TC sensor
15-497-010	Incubator with TC sensor, RH display
15-497-011	Incubator with IR sensor
15-497-012	Incubator with IR sensor, RH display

ACT NOW!

1. Contact your Fisher Scientific sales representative
2. Visit www.thermo.com/protectionpartners to redeem your free accessory

www.thermo.com/protectionpartners



For customer service, call 1-800-766-7000.
To fax an order, use 1-800-926-1166.
To order online: www.fishersci.com

©2008 Thermo Fisher Scientific Inc. All rights reserved. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries.