

## ANALYSIS CERTIFICATE

Manufacturer Name and Address: Thermo Fisher Scientific Oy  
Ratastie 2, FIN-01620 Vantaa, Finland

Product Name : MAGNESIUM

Catalogue no. : 981385

Lot no. : C466

Package : 12 x 20 ml

Storage Conditions : 2...25 °C

Expiry date : 2009-03

### Quality Control

#### Performance test:

<u>Control (lot)</u>	<u>Result (mmol/l)</u>	<u>Target range (mmol /l)</u>
Abtrol (B392)	1,21	1,11 - 1,31
Nortrol (B391)	0,82	0,74 - 0,86

This product with the mentioned lot number was analysed by our authorised quality control department and confirmed that it meets its specifications.

Vantaa 29 / 2 2008



Hannele Laine  
QA Specialist  
Quality Regulatory and Compliance  
Clinical Diagnostics Finland

Konelab 60i (QC)  
SN: A5120360


16.10.2007 14:07 *1kp* *MAGNESIUM 981385 466 2009-03*

Test MG

Accepted - - 17.9.2007 10:41  
Factor 4.714 4.71  
Bias -0.004 -0.005

Coeff. of det. 0.999701 0.998

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors
1	Water	-0.004	0.001	0.000	
2	Water	-0.004	-0.001	0.000	
3	sCal	0.212	1.017	1.030	
4	sCal	0.217	1.043	1.030	<i>LOT B205B</i>