

ANALYSIS CERTIFICATE

Manufacturer Name and Address: Thermo Fisher Scientific Oy
Ratastie 2, FIN-01620 Vantaa, Finland

Product Name : CALCIUM

Catalogue no. : 981367

Lot no. : D341

Package : 6 x 60 ml

Storage Conditions : 2...25 °C

Expiry date : 2010-01

Quality Control

Performance test:

<u>Control (lot)</u>	<u>Result (mmol/l)</u>	<u>Target range (mmol /l)</u>
Abtrol (B392)	2,95	2,87-3,11
Nortrol (B391)	2,11	2,04-2,22

This product with the mentioned lot number was analysed by our authorised quality control department and confirmed that it meets its specifications.

Vantaa 24/7 2008



Hannele Laine
Quality Assurance Specialist
Quality Regulatory and Compliance
Clinical Diagnostics Finland

Konelab 30i (QC)
 SN: P1118753
 Thermo Fisher Scientific (CDX Finland)

23.7.2008 9:55


Test CA 981367 LOT 0341 2010-01

Accepted - - 22.7.2008 7:17

Factor 6.16 6.185
 Bias 0.001 0.001

Coeff. of det. 0.999841 1.000

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors
1	Water	0.001	-0.000	0.000	
2	Water	0.001	0.000	0.000	
3	sCal	0.398	2.448	2.470	
4	sCal	0.405	2.492	2.470	
5	Abtrol (control)	0.474	2.913	2.990	
6	Nortrol (control)	0.346	2.127	2.130	