

ANALYSIS CERTIFICATE

Manufacturer Name and Address: Thermo Fisher Scientific Oy
Ratastie 2, FIN-01620 Vantaa, Finland

Product Name : CALCIUM

Catalogue no. : 981367

Lot no. : E368

Package : 6 x 60 ml

Storage Conditions : 2...25 °C

Expiry date : 2011-05


Quality Control

Performance test:

<u>Control (lot)</u>	<u>Result (mmol/l)</u>	<u>Target range (mmol /l)</u>
Abtrol (D390)	3,18	2,94-3,18
Nortrol (D389)	2,37	2,24-2,48

This product with the mentioned lot number was analysed by our authorised quality control department and confirmed that it meets its specifications.

Vantaa 8 / 10 2009


Hannele Laine
Quality Assurance Specialist
Quality Regulatory and Compliance
Clinical Diagnostics Finland

Konelab 60i (QC)
 SN: A5120360

12.8.2009 10:00 **IKD** **CA E368 2011-04** **981367**


Test CA

Accepted - -

Factor 7.098
 Bias 0.001

Coeff. of det. 0.999984

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors
1	Water	0.001	-0.000	0.000	
2	Water	0.001	0.000	0.000	
3	sCal	0.350	2.477	2.470	B295N
4	sCal	0.348	2.463	2.470	
5	Abtrol (control)	0.444	3.146	3.060	D390A
6	Nortrol (control)	0.322	2.277	2.360	D389A