

## ANALYSIS CERTIFICATE

Manufacturer Name and Address: Thermo Fisher Scientific Oy  
Ratastie 2, FIN-01620 Vantaa, Finland

Product Name : Magnesium

Catalogue no. : 981385

Lot no. : D276

Package : 12 x 20 ml

Storage Conditions : 2...25 °C

Expiry date : 2010-05

**Quality Control  
Performance test:**

<u>Control (lot)</u>	<u>Result (mmol/l)</u>	<u>Target range (mmol/l)</u>
Abtrol (B392)	1,22	1,11 - 1,31
Nortrol (B391)	0,81	0,74 - 0,86

This product with the mentioned lot number was analysed by our authorised quality control department and confirmed that it meets its specifications.

Vantaa 10 / 10 2008



Hannele Laine  
Quality Assurance Specialist  
Quality Regulatory and Compliance  
Clinical Diagnostics Finland

Konelab PRIME 60 KUSTI, SN F4620934□  
 Konelab User

JS

09/06/2008 12:17      981385      D276      2010-05  
 -----


Test      MG

Accepted                      -                      -                      03/06/2008 09:22

Factor                      4.607                      4.62  
 Bias                      -0.004                      -0.003

Coeff. of det.                      0.999740                      1.000

Errors



Calibrator	Response	Calc. con.	Conc.	Errors
1 Water	-0.004	0.001	0.000	B295G
2 Water	-0.004	-0.001	0.000	
3 sCal	0.217	1.018	1.030	
4 sCal	0.222	1.042	1.030	
5 Nortrol (control)	0.172	0.811	0.800	C125D