

ANALYSIS CERTIFICATE

Manufacturer Name and Address: Thermo Fisher Scientific Oy
Ratastie 2, FIN-01620 Vantaa, Finland

Product Name : Magnesium

Catalogue no. : 981385

Lot no. : D113

Package : 12 x 20 ml

Storage Conditions : 2...25 °C


Expiry date : 2010-09

**Quality Control
Performance test:**

<u>Control (lot)</u>	<u>Result (mmol/l)</u>	<u>Target range (mmol/l)</u>
Abtrol (B392)	1,23	1,11 - 1,31
Nortrol (B391)	0,82	0,74 - 0,86

This product with the mentioned lot number was analysed by our authorised quality control department and confirmed that it meets its specifications.

Vantaa 10 / 10 2008


Hannele Laine
Quality Assurance Specialist
Quality Regulatory and Compliance
Clinical Diagnostics Finland

27/02/2008 08:01


/KHe MAGNESIUM, 981385, 0113, 2009-10

Test MG

Accepted - - 12/02/2008 12:38
Factor 4.628 4.65
Bias -0.004 -0.003

Coeff. of det. 0.999624 1.000

Errors



Calibrator	Response	Calc. con.	Conc.	Errors
1 Water	-0.004	0.002	0.000	
2 Water	-0.004	-0.002	0.000	
3 sCal	0.215	1.016	1.030	
4 sCal	0.221	1.044	1.030	<i>0.295</i>