



Advantage – advanced analysis software for Thermo Scientific DCA RADIANT




Thermo Scientific Advantage is a stand-alone Windows software package designed for use with the DCA RADIANT 300 series, (Dynamic Contact Angle) systems operated under the WinDCA software. Advantage provides advance analysis of DCA RADIANT files allowing the calculation of the solid surface energy, the liquid surface tension components, the wicking rate, and the wettability of powder samples.

- Opens standard WinDCA format files
- Quick multiple file selection and loading
- Standard Microsoft Office interface
- Flexible and easy-to-use
- Runs under Windows 2000




Analysis modules

- Zisman plot
- Geometric mean
- Harmonic mean
- Lewis Acid/Base
- Liquid surface tension components
- Wicking rate
- Powder analysis



Graphics

- Advanced graphing capabilities
- Easy to use interface
- Mouse operated roll and zoom



Results

- Complete calculations
- Easy file management
- Printing of reports
- Windows copy and paste functions

**Thermo Fisher Scientific
Process Instruments**

International/Germany
Dieselstr. 4,
76227 Karlsruhe
Tel. +49(0)721 40 94-444
info.mc.de@thermofisher.com

Benelux
Tel. +31 (0) 76 5 87 98 88
info.mc.nl@thermofisher.com

China
Tel. +86 (21) 68 65 45 88
info.mc.china@thermofisher.com

France
Tel. +33 (0) 1 60 92 48 00
info.mc.fr@thermofisher.com

India
Tel. +91 (22) 27 78 11 06
info.mc.in@thermofisher.com

United Kingdom
Tel. +44 (0) 1785 81 36 48
info.mc.uk@thermofisher.com

USA
Tel. 603 436 9444
info.mc.us@thermofisher.com

www.thermo.com/mc

C002_13.08.2008

© 2007/01 Thermo Fisher Scientific.
All rights reserved. This document is
for informational purposes only and is
subject to change without notice.