



Selecting the Correct Phase by Industry Application

Industry	Application	Recommended Sorbent
Forensics	Clean up of acidic drugs from biological fluids	HyperSep Verify-AX
	Clean up of basic drugs from biological fluids	HyperSep Verify-CX
Pharmaceutical	Clean up of polar metabolites from plasma / serum	HyperSep Hypercarb
	Clean up of basic drugs	HyperSep Verify-CX
	Clean up of barbiturates from biological fluids	HyperSep C18
	Clean up of antineoplastic agents from biological fluids	HyperSep C18
	Clean up of beta blockers from biological fluids	HyperSep C18, HyperSep CN
	Clean up of benzodiazepines from biological fluids	HyperSep C18
	Clean up of Diuretics from biological fluids	HyperSep C18
	Environmental	Clean up of PAH's from water
Clean up of PAH's from oil		HyperSep FL
Clean up of PAH & PCB's from water		HyperSep C18
Clean up of pesticides from water		HyperSep C18
Clean up of aromatic hydrocarbons from water		HyperSep C18
Clean up of trace metal elements from water		HyperSep PH, HyperSep C18, HyperSep NH2
Biochemical	Clean up of catecholamines from biological fluids	HyperSep SAX, HyperSep SCX
	Clean up of amino acid from solution	HyperSep SCX
	Clean up of preservatives from cosmetics	HyperSep C18, HyperSep SCX
	Clean up of mycotoxins from foodstuff	HyperSep SI, HyperSep C18, HyperSep PH
	Clean up of nucleosides from solution	HyperSep C18
	Clean up of steroids from biological fluids	HyperSep C18

