

ANALYSIS CERTIFICATE

Manufacturer Name and Address: Thermo Fisher Scientific Oy
Ratastie 2, FIN-01620 Vantaa, Finland

Product Name : CALCIUM

Catalogue no. : 981367

Lot no. : D497

Package : 6 x 60 ml

Storage Conditions : 2...25 °C

Expiry date : 2010-01


Quality Control

Performance test:

<u>Control (lot)</u>	<u>Result (mmol/l)</u>	<u>Target range (mmol /l)</u>
Abtrol (B392)	2,95	2,87-3,11
Nortrol (C125)	2,09	2,04-2,22

This product with the mentioned lot number was analysed by our authorised quality control department and confirmed that it meets its specifications.

Vantaa 7 / 11 2008



Hannele Laine
Quality Assurance Specialist
Quality Regulatory and Compliance
Clinical Diagnostics Finland

29/10/2008 10:23

381367 0497 2010-01


Test CA

Accepted - -

Factor 6.149
 Bias 0.001

Coeff. of det. 0.999956

Errors



Calibrator	Response	Calc. con.	Conc.	Errors
1 Water	0.001	0.002	0.000	
2 Water	0.001	-0.002	0.000	
3 sCal	0.405	2.481	2.470	
4 sCal	0.401	2.459	2.470	
5 Abtrol (control)	0.471	2.888	2.990	
6 Nortrol (control)	0.344	2.110	2.130	