

## ANALYSIS CERTIFICATE

Manufacturer Name and Address: Thermo Fisher Scientific Oy  
Ratastie 2, FIN-01620 Vantaa, Finland

Product Name : Magnesium

Catalogue no. : 981385

Lot no. : D495

Package : 12 x 20 ml

Storage Conditions : 2...25 °C


Expiry date : 2010-05

**Quality Control  
Performance test:**

<u>Control (lot)</u>	<u>Result (mmol/l)</u>	<u>Target range (mmol/l)</u>
Abtrol (B392) dil1+3	1,20	1,11 - 1,31
Nortrol (B391)	0,82	0,74 - 0,86

This product with the mentioned lot number was analysed by our authorised quality control department and confirmed that it meets its specifications.

Vantaa 7 / 11 2008

  
Hannele Laine  
Quality Assurance Specialist  
Quality Regulatory and Compliance  
Clinical Diagnostics Finland


27/10/2008 11:58 / *KMC* *981385* *LOT 0495 exp 2010-05*  
 -----

Test MG

Accepted - -  
 Factor 4.544  
 Bias -0.006

Coeff. of det. 0.999716

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors
1	Water	-0.006	0.001	0.000	
2	Water	-0.006	-0.001	0.000	
3	sCal	0.224	1.042	1.030	
4	sCal	0.218	1.018	1.030	
5	Nortrol (control)	0.173	0.814	0.800	