

Material Safety Data Sheet

Phoenix Blue Nuclear Stain-Solution B

Section 1 - Chemical Product and Company Identification

MSDS Name:

Phoenix Blue Nuclear Stain - Solution B

Catalog Numbers:

7214

Synonyms:

None Known.

Company Identification:
 Richard Allan Scientific
 4481 Campus Drive
 Kalamazoo, MI 49008
Company Phone Number:

800-522-7270

Emergency Phone Number:

800-424-9300

CHEMTREC Phone Number, US:

(800) 424-9300

CHEMTREC Phone Number, Europe:

(202) 483-7616

Section 2 - Composition, Information on Ingredients

CAS#	Chemical Name:	Percent	EINECS/ ELINCS	Hazard Symbols	Risk Phrases
7732-18-5	Deionized Water	95-100	231-791-2		
Not available	Proclin 300	<1.0	Not available		
	Proprietary Ingredient	<1.0	216-346-2		

Section 3 - Hazards Identification

EMERGENCY OVERVIEW

Appearance: Blue liquid

Caution! May cause eye, skin, and respiratory tract irritation. The toxicological properties of this material have not been fully investigated.

Target Organs: None known

Potential Health Effects**Eye:**

May cause eye irritation.

Material Safety Data Sheet

Phoenix Blue Nuclear Stain-Solution B

Skin:

May cause skin irritation. May be harmful if absorbed through the skin.

Ingestion:

May cause gastrointestinal irritation with nausea, vomiting and diarrhea. May be harmful if swallowed.

Inhalation:

May cause respiratory tract irritation. May be harmful if inhaled.

Chronic:

No information found.

Section 4 - First Aid Measures

Eyes:

Immediately flush eyes with plenty of water for at least 15 minutes, occasionally lifting the upper and lower eyelids. Get medical aid.

Skin:

Get medical aid. Flush skin with plenty of water for at least 15 minutes while removing contaminated clothing and shoes.

Ingestion:

Get medical aid. Do NOT induce vomiting. If conscious and alert, rinse mouth and drink 2-4 cupfuls of milk or water.

Inhalation:

Remove from exposure and move to fresh air immediately. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get medical aid.

Notes to Physician:

Treat symptomatically and supportively.

Section 5 - Fire Fighting Measures

General Information:

As in any fire, wear a self-contained breathing apparatus in pressure-demand, MSHA/NIOSH (approved or equivalent), and full protective gear. During a fire, irritating and highly toxic gases may be generated by thermal decomposition or combustion.

Extinguishing Media:

Use water spray, dry chemical, carbon dioxide, or chemical foam.

Autoignition Temperature:

Not available

Explosion Limits:

Lower: Not available Upper: Not available

Flash Point:

Not available

NFPA Rating:

(estimated) Health: 1; Flammability: 0; Instability: 0

Material Safety Data Sheet

Phoenix Blue Nuclear Stain-Solution B

Section 6 - Accidental Release Measures

General Information:

Use proper personal protective equipment as indicated in Section 8.

Spills/Leaks:

Absorb spill with inert material (e.g. vermiculite, sand or earth), then place in suitable container. Clean up spills immediately, observing precautions in the Protective Equipment section. Provide ventilation.

Section 7 - Handling and Storage

Handling:

Avoid breathing dust, mist, or vapor. Avoid contact with eyes, skin, and clothing. Keep container tightly closed. Avoid ingestion and inhalation. Use with adequate ventilation.

Storage:

Store in a cool, dry place. Store in a tightly closed container.

Section 8 - Exposure Controls, Personal Protection

Engineering Controls:

Facilities storing or utilizing this material should be equipped with an eyewash facility and a safety shower. Use adequate ventilation to keep airborne concentrations low.

Exposure Limits

Chemical Name:	ACGIH	NIOSH	OSHA
Deionized Water	None listed	None listed	None listed
Proclin 300	None listed	None listed	None listed
Proprietary Ingredient	None listed	None listed	None listed

OSHA Vacated PELs
Personal Protective Equipment
Eyes:

Wear appropriate protective eyeglasses or chemical safety goggles as described by OSHA's eye and face protection regulations in 29 CFR 1910.133 or European Standard EN166.

Skin:

Wear appropriate protective gloves to prevent skin exposure.

Clothing:

Wear appropriate protective clothing to prevent skin exposure.

Respirators:

Follow the OSHA respirator regulations found in 29 CFR 1910.134 or European Standard EN 149. Use a NIOSH/MSHA or European Standard EN 149 approved respirator if exposure limits are exceeded or if irritation or other symptoms are experienced.

Material Safety Data Sheet
Phoenix Blue Nuclear Stain-Solution B**Section 9 - Physical and Chemical Properties**

Physical State: Liquid
Color: Blue
Odor: No information found
pH: No information found
Vapor Pressure: No information found
Vapor Density: No information found
Evaporation Rate: No information found
Viscosity: No information found
Boiling Point: No information found
Freezing/Melting Point: No information found
Decomposition Temperature: No information found
Solubility in water: No information found
Specific Gravity/Density: No information found
Molecular Formula: Solution
Molecular Weight: No information found

Section 10 - Stability and Reactivity**Chemical Stability:**

Stable under normal temperatures and pressures.

Conditions to Avoid:

Incompatible materials, excess heat

Incompatibilities with Other Materials

Strong oxidizing agents.

Hazardous Decomposition Products

Carbon monoxide, carbon dioxide

Hazardous Polymerization

Will not occur.

Section 11 - Toxicological Information**RTECS:**

CAS# 7732-18-5: ZC0110000
Proprietary Ingredient unlisted

LD50/LC50:

CAS# 7732-18-5:
Oral, rat: LD50 = >90 mL/kg.
Proprietary Ingredient:
No information found

Carcinogenicity:

CAS# 7732-18-5: Not listed as a carcinogen by ACGIH, IARC, NTP, or CA Prop 65.

Material Safety Data Sheet
Phoenix Blue Nuclear Stain-Solution B

Proprietary Ingredient: Not listed as a carcinogen by ACGIH, IARC, NTP, or CA Prop 65.

Epidemiology:

No information found

Teratogenicity:

No information found

Reproductive:

No information found

Mutagenicity:

No information found

Neurotoxicity:

No information found

Other:

The toxicological properties have not been fully investigated.

Section 12 - Ecological Information

No information found

Section 13 - Disposal Considerations

Chemical waste generators must determine whether a discarded chemical is classified as a hazardous waste. US EPA guidelines for the classification determination are listed in 40 CFR Part 261.3. Additionally, waste generators must consult state and local hazardous waste regulations to ensure complete and accurate classification.

RCRA P Series Wastes

None of the components are on this list.

RCRA U Series Wastes

None of the components are on this list.

Section 14 - Transport Information**US DOT****Canadian TDG**

Proper Shipping Name: Not Regulated

Not Regulated

Hazard Class:

UN Number:

Packing Group:

Material Safety Data Sheet

Phoenix Blue Nuclear Stain-Solution B

State Right to Know

No information found

California Prop 65

None of the components are on this list.

California No Significant Risk Level

None of the components are on this list.

None of the components are on this list.

None of the components are on this list.

European/International Regulations**European Labelling in Accordance with EC Directives:**

Hazard Symbols: None listed

Risk Phrases: None listed

Safety Phrases: S 24/25 Avoid contact with skin and eyes.

WGK (Water Danger/Protection)

No information found

United Kingdom Occupational Exposure Limits

No information found

United Kingdom Maximum Exposure Limits

No information found

Canadian DSL/NDSL

CAS# 7732-18-5 is listed on Canada's DSL List.

Canadian WHMIS Classifications

No information found

This product has been classified in accordance with the hazard criteria of the Controlled Products Regulations and the MSDS contains all of the information required by those regulations.

Canadian Ingredient Disclosure List

No information found

Section 16 - Other Information

No information found

MSDS Creation Date: April 15, 2008

Revision Date: April 15, 2008

Revisions were made in Sections:

6, 9, 14, 16

This MSDS is intended for review and guidance in the receipt, storage, handling, use and disposal of product purchased from us, and for no other purpose. Use this product only as directed and in accordance with applicable instructions and warnings provided with the product. Please consult your institution's policies regarding use of this product. If you have obtained this MSDS other than in connection with the supply of this product from us, this MSDS should be consulted for general information only, and should not be relied upon for any purpose. As with the use of all hazardous materials, you should in all instances follow the guidance of the MSDS provided or available with the specific product purchased.